

Using EpiPen®/EpiPen® Jr is as easy as 1-2-3

1. Remove yellow or green cap from carrying case
 - Grasp unit with black tip pointing downward
 - Pull off grey activator cap



2. Jab black tip firmly into outer thigh so it "clicks" AND HOLD on thigh approximately 10 seconds
 - Massage injected area for 10 seconds



3. Seek medical attention



After Using EpiPen®/EpiPen® Jr Follow 3 Easy Safety Steps:

1. Carefully place used auto-injector, needle-end first, into storage tube



2. Screw cap of carrying case on completely.
 - This automatically bends needle back and secures pen so it won't fall out of tube



3. Take unit with you to hospital Emergency Department



HOW TO ADMINISTER A TWINJECT

STEP 1 MAKE SURE THAT THE MEDICINE IS READY

Look at Twinject™ 0.3 mg or Twinject™ 0.15 mg regularly. It may not work if medicine looks cloudy (has particles), pinkish, or more than slightly yellow, or if the expiration date has passed.

In the event of a life-threatening allergic reaction, you should use an out of date product, if that is all you have. Do NOT remove the GREEN or the RED cap until you are ready to use.

STEP 2 FIRST DOSE

1. PULL off GREEN end cap to see a GREY cap. **Never put thumb, finger, or hand over the GREY cap.**
2. PULL off RED end cap.
3. Place GREY cap against mid-thigh (can go through clothes).
4. Press down firmly until auto-injector activates-hold while slowly counting to ten.
5. Remove auto-injector and check the GREY cap; if needle is exposed, you received the dose. If not, repeat #3 and #4 under step 2. Prepare for second dose.



Pull off green end cap



Press down firmly until auto-injector activates

Get emergency medical help right away.

STEP 3 PREPARE FOR A SECOND DOSE

1. Unscrew and remove GREY cap. **Beware of exposed needle.**
2. Holding **BLUE** hub at needle base, remove syringe from barrel.
3. Slide **YELLOW** (Twinject 0.3 mg) or **ORANGE** (Twinject 0.15 mg) collar off plunger.
4. **PAUSE HERE.** If symptoms have not improved in approximately 10 minutes since first injection, proceed with Step 4.



Holding blue hub...



Slide yellow collar

STEP 4 INJECT SECOND DOSE

1. Insert needle into mid-thigh.
2. Push plunger down completely.

Get emergency medical help right away.



Insert needle into mid-thigh

SYRINGE DISPOSAL: Re-insert syringe, needle first, into blue case. Return it to your physician or pharmacist for proper disposal. **Do not throw away in trash.**